**FISBA** Innovators  
in Photonics

# Technical Instruction

**FISBA READYBeam™**  
**Compact multi-wavelength laser source**

January 2020

# **Content**

<b>1</b>	<b>Safety</b>	<b>4</b>
1.1	Explanation of symbols	4
1.1.1	General safety instructions	4
1.1.2	Laser safety	4
1.1.3	Laser area	4
1.2	Sécurité	6
1.2.1	Explication des symboles	6
1.2.2	Consignes de sécurité générales	6
1.2.3	Sécurité du laser	6
1.2.4	Zone laser	6
1.3	Warning and information labels	8
1.4	Correct and specified use of the unit	8
1.5	Warranty conditions	9
1.5.1	Limits of the warranty	9
<b>2</b>	<b>Technical system</b>	<b>9</b>
2.1	Laser module	9
2.2	Installation and commissioning	10
2.2.1	Unpacking and checking for visible damage	10
2.2.2	Installation	10
2.2.3	Operating the FISBA READYBeam™	10
2.2.4	Digital operating with the software application	10
2.2.5	Software interface	11
2.2.6	Analog controlled operation	11

# **1 Safety**

## **1.1 Explanation of symbols**

This symbol is used for all points in these operating instructions to which special attention is to be paid so that all directives, regulations, instructions and the correct work sequence are observed, and to prevent minor or serious damage to the laser system or plant.

### **1.1.1 General safety instructions**

The operating instructions and the safety instructions are to be read and observed prior to installation and commissioning! The operating instructions must be read, understood and followed by the operators responsible for the unit. Along with the instructions in these operating instructions, also observe the generally applicable safety and accident prevention regulations! All work on installation and commissioning as well as all maintenance is to be performed by appropriately qualified personnel. National accident prevention regulations, EN and IEC standards are to be observed. Appropriately qualified personnel in the context of these basic safety instructions are persons who are familiar with the fitting, installation, placing in operation and operation of the product and have qualifications appropriate to their task.

### **1.1.2 Laser safety**

The unit is a class 3B laser device. Both the direct beam and its reflections from diffuse reflective surfaces are dangerous. The unit emits strong power beams in the visible spectral range (approx. 400 nm to 700 nm, for more detailed information see the data sheet, page 15). The radiation can cause irreversible damage to the eyes if the necessary protective measures are not taken.

### **1.1.3 Laser area**

The laser area is the area in which the values for the maximum permissible irradiation can be exceeded. Here the possibility of unintentional deflection of the laser beam is also to be taken into account. If operational equipment other than that given in this manual or another procedure is used, dangerous exposure to the radiation may result.

## **Laser system screened inside a protective cover**

By using suitable screening it can be ensured that people are not exposed to the laser radiation. In this way the laser class is reduced and it is possible to work with the closed arrangement without additional limitations or further measures.

However, in this case the following requirements apply to the protective cover:

- The protective cover must be appropriate to adequately shield the laser radiation. It is therefore necessary to use laser safety glass or laser safety film for viewing windows.
- The cover itself as well as any doors or flaps in the cover must be protected using the interlock circuit that automatically shuts down and inhibits the laser immediately on intentional or unintentional opening.

If the interlock circuit is disabled for service or maintenance work on the system with the cover open, the laser area is enlarged again and the regulations in the next section «Presence of persons in the laser area» apply.

## **Presence of persons in the laser area**

In the laser area, people are subject to hazards, in particular ocular. Therefore the following measures must be taken and the following safety regulations observed:

- The room must have an emergency exit.
- The room must be optically screened from the environment (e.g. laser safety film on the windows).
- The laser area is to be kept as small as possible, to be bounded by suitable screening and protected against access by unauthorized persons. The number of persons in the laser area should be reduced to the minimum.
- There must be a warning device on all entrances to the room that indicates the laser radiation hazard.
- If the entry doors are not protected using the interlock circuits, opening the doors easily from the outside must be impossible to prevent thoughtless entry.
- All persons who are in the laser area during operation must be informed about the dangers of the laser radiation.
- Laser safety glasses/goggles must be worn in the laser area. The operating organisation must ensure that suitable safety glasses/goggles are available. The glasses/goggles are selected in accordance with DIN EN 207. Please see the data sheet on the unit for the actual wavelength of your laser. Hereinafter, the term safety glasses is always used to refer to suitable laser safety glasses/goggles. **Caution:** in some circumstances your laser safety glasses/goggles may not protect you from a powerful red pilot beam.

- There must not be any potentially explosive substances in the laser area. Easily inflammable substances may catch fire.
- Glossy reflective equipment must be either removed from the laser area or covered. Windows and reflecting walls are also to be covered with material with low inflammability.

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### 1.3 Warning and information labels

The following warning labels are fitted to laser and processing head:

<b>VISIBLE LASER RADIATION AVOID EXPOSURE TO BEAM</b>	<b>VISIBLE LASER RADIATION AVOID EXPOSURE TO BEAM</b>
<b>CLASS 3B LASER PRODUCT</b>	<b>CLASS 3B LASER PRODUCT</b>
$\lambda[\text{nm}]$	max. power [nW]
400–410	210
486–490	55
632–643	155
EN 60825-1:2014	EN 60825-1:2014
 	<b>DANGER</b>
This product does not confirm to 21CFR 1040 and EN 60825-1, it is intended to be integrated in OEM equipment only.	Laser radiation. Avoid direct exposure to Beam
	 400–700 nm      Class 111b Max. 300 mW      laser product.

### 1.4 Correct and specified use of the unit

It is forbidden to commission the laser system until it has been ensured that the machine or plant in which the laser system has been fitted complies with the

applicable safety regulations. The organisation operating the plant is solely liable for any damage caused by incorrect use of the unit. The organisation operating the system is responsible for the correct conditions in relation to operation, maintenance and repair. The only personnel to be tasked with operation, maintenance and repair are personnel who have been adequately instructed in the function of the laser system and who have also been instructed specifically on the dangers and risks of incorrect handling of the system.

## 1.5    **Warranty conditions**

The warranty period is 12 months from delivery. The warranty covers the entire laser module. The warranty will become void on the

- unauthorized opening of the unit's components
- operation of the unit in an unauthorized configuration
- improper use, storage or transport (e.g. vibration, temperature shock, the action of frost)

### 1.5.1    **Limits of the warranty**

No warranty of the suitability of the product for specific applications is provided. FISBA is not liable for indirect, direct or consequential damages caused by the use of this product.


## 2    **Technical system**

The laser module emits laser radiation with different wavelengths. Combination of different laser diodes can emit radiation in the UV, VIS and/or NIR range. The correct combination is documented within the data sheet, see page 15.

## 2.1 Laser module

Dimensions of the laser module in mm:

Connector: Samtec SFSD-10-28G24.00SR

Pin assignment: documented in the appendix (Table 1 on page 13)

## 2.2 Installation and commissioning

Note that the laser module with electronic driver is not a product for direct use. Security elements, such as interlock, and emergency switch are not included in the laser module. If the laser module is used as part of system or product, then the responsibility for safety lies with the product development team of the product.

When operating the laser module, the laser safety regulations must be observed.

### 2.2.1 Unpacking and checking for visible damage

Check the completeness of the delivery and check that all the items supplied are in good condition. In case of errors please contact FISBA AG ([readybeam@fisba.com](mailto:readybeam@fisba.com)). If possible, keep the packaging in which the laser system was supplied. It will then be possible to pack the unit in its original packaging and transport it safely in case of repair.

## **2.2.2 Installation**

The laser module can be mounted with screws on a flat bottom plate. The bottom plate should be cooled and the heat transfer should be at least 7W to hold the module on a stable operating temperature. The laser module READYBeam is TEC-controlled.

Incorrect cooling effects the laser power and can destroy the laser diodes. During assembly of the laser module on a plane cooling surface, attention is required so that no tension can apply to the casing. Tension influences the beam quality very strongly. The connection must be handled carefully.


## **2.2.3 Operating the FISBA READYBeam™**

- The safety regulations must be carefully observed
- A power supply between +12 and +24V with max. 4 Ampere is needed to start the READYBeam
- FISBA READYBeam™ can be started digital with the software application and RS485 interface or analog based controlled

## **2.2.4 Digital operating with the software application**

- Visit [fisba.com/readybeam-software](http://fisba.com/readybeam-software) to download the software
- Unpack the files from downloaded folder
- Start «RGB Service Software Setup.msi» and follow the installation instructions
- You need an USB - RS485 interface cable for the digital operating mode  
See the details in the appendix (Figure 1 on page 14)
- Connect the USB plug to your laptop and the RS485 plug to the laser module interface cable
- Connect the power supply to the interface cable
- Connect the interface cable to the laser module connector
- Start the READYBeam software
- Now it is possible to control the lasers via the software

## 2.2.5 Software interface

1. Set **Modulation Mode**: choose «Digital»
2. Set **Power Value 1** (red), **Power Value 2** (green) and **Power Value 3** (blue): values in % of the max. laser power
3. Switch the Laser on/off with **Digital Enable 1, 2 or 3**
4. Enter the command with «**Write Config**» Button

The READYBeam software is provided for download at [fisba.com/readybeam-software](http://fisba.com/readybeam-software)

## 2.2.6 Analog controlled operation

- Connect the power supply and signals according the pin configuration on Table 1, page 13
- Connect the Interface cable at READYBeam
- Switch ON +12/+24V Power

## Appendix

Pin configuration and power values

Connector: Samtec SFSD-10-28G24.00SR  
Enable Laser: low = 0V, high = 3.3V  
Analog In: 0V ... 3.3V → 10% ... 100% of Laser power  
Power: +12V ... +24V max. 4A

Pin 1 is marked on the connector with arrow (not notch)	
Pin 1	Analog In- Laser red
Pin 2	Analog In+ Laser red
Pin 3	Analog In- Laser green
Pin 4	Analog In+ Laser green
Pin 5	Analog In- Laser blue
Pin 6	Analog In+ Laser blue
Pin 7	Gnd
Pin 8	Enable Laser red
Pin 9	Enable Laser green
Pin 10	Enable Laser blue
Pin 11	RS 485 B
Pin 12	RS 485 A
Pin 13	Gnd
Pin 14	Gnd
Pin 15	Gnd
Pin 16	Gnd
Pin 17	Power +12V ... +24V
Pin 18	Power +12V ... +24V
Pin 19	Power +12V ... +24V
Pin 20	Power +12V ... +24V

Table 1

Figure 1

## Data sheet FISBA READYBeam™

### Technical Specifications

Optical parameters	bio			ind					
Wavelengths	405 nm	488 nm	638 nm	450 nm	520 nm	660 nm			
Min. optical power from fiber	≥ 40 mW	≥ 30 mW	≥ 40 mW	≥ 40 mW	≥ 30 mW	≥ 40 mW			
NA fiber	0.12								
Laser class	3B								
Fiber connector	FC/APC								
Fiber type	Single-mode with or without polarization maintaining properties								
Electrical parameters									
Power supply	12 V ... 24 V								
Modulation	ON / OFF Enable up to 1 MHz Set power input bandwidth up to 20 kHz (0..3.3 V -> 10%..100%)								
Operation mode	cw / modulated								
Digital interface	RS 485 (ON / OFF Modulation up to 100 Hz)								
Mechanical parameters									
Dimensions length x width x height	77 mm x 40 mm x 37 mm								
Weight	≤ 200 g								
Environmental conditions									
Operation temperature	15°C – 40°C								
Storage temperature	-10°C ... +60°C								
Thermal properties	TEC controlled								
Heat dissipation	~ 7 watts								

<b>FISBA AG</b>	Rorschacher Str. 268 9016 St.Gallen Switzerland
<b>FISBA Photonics GmbH</b>	Schwarzschildstrasse 10 12489 Berlin Germany
<b>FISBA LLC</b>	6296 E. Grant Rd Suite 150 Tucson, AZ 85712 United States
<b>FISBA (Shanghai) Co., Ltd.</b>	9F-4 World Square Building, No. 855 Pudong South Rd. Shanghai China